


HONKON 808CUTE

808nm Diode Laser Portable
Permanent Hair Removal Machine

**Premium German Technology
Hair Reduction Laser**

800W Micro-channel diode laser
Big spot size: 12mm*23mm
6mm treatment tip for small size
Portable SHR permanent hair removal machine
Safe/Comfortable/Fast/Painless

**APPLICATIONS:
Hair Removal
Skin Rejuvenation**



FDA | CE



+91-94460-01472
+91-495-2742210



info@evolutionmed.com
www.evolutionmed.com

Evolution-Med is the authorised distributors of Beijing HONKON Technologies Co. Ltd. in Kerala and Tamil Nadu. HONKON is a leading professional Medical & Beauty Equipment manufacturer engaged in beauty industry more than 21 years and started operations in India in 2001.

HONKON Service Centers: MUMBAI, DELHI, BANGALURU, HYDRABAD, CHENNAI & COCHIN



HONKON

808CUTE


808nm Diode Laser
Portable Permanent
Hair Removal Machine




Treatment Principle



The laser selectively acts on the melanin in the hair follicle, which destroys the germinal region in the hair warm.



Natural hair shedding, to achieve the purpose of hair removal



Stimulate collagen regeneration, reduce the pores make the skin tight smooth at the same time.

Features & Advantages

- Fashion type & Cost-effective, convenient for transportation.
- Professional and effective on all skin colors hair removal treatment.
- 800w 808 cute achieves the safest and most comprehensive hair removal.
- Micro channel for cooling system is more comfortable and safe.
- 12.1 inches latest touch display screen big and convenient to operate.
- Strong cooling systems, water cooling+ air cooling+ Semiconductor Cooling.
- Big Spot Size: 12x23mm and 1-10 HZ high frequency in motion treatment can be realized to make sure the machine can do large area hair removal treatment in short time.



+91-94460-01472

+91-495-2742210

info@evolutionmed.com

www.evolutionmed.com

UNIQUE ADVANTAGES

Laser power upto 800W

808CUTE diode laser beauty machine, 800W highpower, narrow pulse width, ultra-ice hair removal, high speed and accurate.

Large Spot and Deep Penetration

808CUTE Large spot size (12mm x 23mm) Deep laser penetration ensures more secure hair removal effect.

Adopting advanced micro channel technology

The most cost-effective portable hair removal equipment with micro channel technology worldwide realizes diode laser freezing painless hair removal.

Portable - Convenient to put it in clinic & carry

Tinny and portable diode laser hair removal machine with powerful micro channel technology is more suitable for transportation and carrying for mobile medical treatments.

BEFORE & AFTER

