




CleanN

N Times Cleaning Power
5 Steps to Nourish One's Skin



CleanN

N Times Cleaning Power
5 Steps to Nourish One's Skin

CleanN is a non-invasive cleansing and nourishing advanced customized beauty equipment with a full set of popular solutions. CleanN takes the new generation of dual-power fluid vortex technology as the core. With the aid of vacuum negative pressure and three-zone vortex probe, CleanN works with natural botanical extracts, and can “inhale” and “conduct” without damaging the skin, resulting in N times the cleansing power to purify, moisturize and rejuvenate the skin. CleanN also integrates five technologies perfectly, including new vibration superconductivity, sonic spade, hot and cold superconductivity, water oxygen cleansing and high-pressure water replenishment. It combines with a wide range of products in a non-invasive way to meet the advanced customized needs of different beauty seekers for cleansing and nourishing their faces and skin.



Appearance & Design

Can't let go of
the black and
gold design.

**Champagne black gold
ultra-thin bezel design.**

High-end WYSIWYG (what you
see is what you get)

**15-inch high-definition
visualization touch
integrated screen.**

More intuitive visualization and
lighter operation.

**Industrial design
appearance.**

Ultra-smooth lines
Ultra-metallic materials
Ultra-refined craftsmanship

**16770000 color
adjustable ambient
lighting**

Colorful atmosphere
starts right now



Full ergonomic
design for comfort
and safety.

**Aluminum handle
design.**

Ergonomic Smoother and
lighter operation.

**Full-size aluminum
handrail.**

Pushing and pulling more
smoothly.

**Plug-and-play bottle
design.**

Saving time and effort.

**Integrated reinforced
glass.**

Texture within reach.

**Large-sized silent
wheels.**


More Silence in experience.

5 Steps to Nourish One's Skin

Cleansing and nourishing, once every two weeks



1 Promoting Absorption




2 Exfoliating



3 Soothing & Improving



4 Cleansing & Revitalizing



5 Hydrating & Rejuvenating

1. Promoting Absorption

Skin Rejuvenation & Moisturizing

Vibrating Superconducting Handles

Enhancing absorption of products

1.1MHz/s vibration frequency

2. Smooth & comfortable exfoliation

Sonic Spade Handles

30KHz Ultra-high frequency vibration

High-Frequency Vibration Sonic Technology

3. Dual Temperature to Sooth & Improve

36-46°C High-temperature

0-6°C Low-temperature

Hot and Cold Temperature Superconducting Technology

Hot and Cold Superconducting Handles

4. Cleansing & Revitalizing

New Chemical Butterfly Water Oxygen Technology

230m/s Ultra-high-speed spraying

50-80um Ultra-small diameter particles

Water Oxygenation Cleaning Handles

Epidermis

Dermis

Hypodermis Tissue

5. Hydrating & Rejuvenating

Non-invasive High-pressure Jet Hydration

50-850KPA Precise pressure control

450m/s High-speed injection

Nano-sized molecules after product automatization

Cuticle

Spinose

Basal Layer

High Pressure Hydration Handle

N Times Cleaning Power


The new generation of dual power fluid vortex technology



With the aid of vacuum negative pressure and three-zone vortex probe, the new generation of dual-power fluid vortex technology can fit the skin accurately from different angles, suck the dirt in the pores efficiently without damaging the skin and introduce the natural plant extracts into the skin at the same time. It combines "inhalation" and "conduction", so as to achieve cleansing and nourishing effect to purify, moisturize and rejuvenate the skin.

Three-Zone Vortex Probe

Exfoliating Soft Blade
Dirt Inhalation
Essence Conduction



The three-zone vortex probe adopts the three-zone design of "exfoliating soft blade, dirt inhalation and essence conduction". At the same time, with the rounded and non-sharp process, it can fit the skin accurately at different angles, and does not irritate or damage the skin while using.

N natural plant extracts 3 steps to cleanse and nourish one's skin

Non-intrusive combination of "inhalation" and "conduction" with natural plant extracts. N times the cleansing power to purify, moisturize and rejuvenate the skin.

1 Purifying & Revitalizing

Clarifying Nib + Clarifying Serum

Cleansing the epidermis gently and opening the skin's channels to make the skin clearer and more energized.

2 Hydrating & Moisturizing

Hydrating Nib + Hydrating Serum

Absorbing dirty waste deeply from pores, regulating oil balance and keeping skin radiant and moisturized.

3 Resurfacing & Rejuvenating

Glow Nib + Glow Serum

Repairing serum effectively penetrates deep into the skin to enhance the skin's protective power and care for facial hydration.



+91-94460-01470
+91-94460-01472



info@evolutionmed.com
www.evolutionmed.com