

SHE-MSP001

Magshape

Multi functional
Body Shaping
& Muscle Building
Instrument



Two
Technology
in One
**HI-EMT
& RF**

30 Mins Treatment = **20.000** Workouts

Combined 30 Min. Treatment

- Less 30%
FAT

More 25%
MUSCLE **+**




info@evolutionmed.com
www.evolutionmed.com

Magshape

Multi Functional Body Shaping & Muscle Building Instrument

Treatment Principle



Through HI-EMT (high energy focused electromagnetic wave) technology, autologous muscles can realize the ultimate training of continuous expansion and contraction, so as to deeply reshape the internal structure of muscles. It promotes the growth of myofibrils (muscle enlargement) to produce new protein chains and muscle fibers (muscle hyperplasia) and increases muscle density and volume through training.

Under Hi-EMT technology, the 100% limit muscle contraction movement can cause a large amount of fatty acid decomposition, fatty acid from the triglyceride decomposition outflow, in the fat cell a large amount of accumulation. Fatty acid concentration is too high, allowing adipocyte apoptosis, in a few weeks by the body normal metabolism out of the body.

So HI-EMT can strengthen and increase muscle, at the same time accomplish reducing fat effectively.

A specific range of frequencies that do not allow the muscle to relax between successive stimuli, forcing the muscle to remain contractive for several seconds. Muscle formation occurs when muscle tissue is repeatedly exposed to these high load conditions.

Two Technology in One HI-EMT & RF

- HIEMT Pro is world-first treatment that combines HIEMT technology with RF.
- HIEMT: Electromagnetic fields trigger thousands of muscle contractions with high intensity, which forces the muscle fibres and cells to multiply and grow at a much faster rate, helping patients to significantly increase their muscle volume.
- RF: While radiofrequency heat can burn and destroy fat cells. The radiofrequency heating prepares muscles for exposure to stress, similar to what a warm up activity does before any workout. It boosts blood supply to muscles so they build up stronger and faster.
- It also heats layers of subcutaneous fat to a level that causes apoptosis and physically break down. Then damaged fat cells will be slowly removed by body metabolism in a natural way, leading to permanent fat removal.
- The radiofrequency process also helps to improve skin laxity while tightening up the skin at the same time.



Features & Advantages

- Building muscles and burn fat together.
- Non-invasive hip lifting.
- Suitable for everyone-no anesthesia-no surgery.
- Short-term treatment, one time 30 minutes.
- Only 4 sessions are required between 5-10 days.
- Feels like an intensive exercise.
- Safe without downtime.
- The result is immediate, but it will be noticeable after two to four weeks.
- Muscle mass increased by an average of 25%.
- Average fat loss 30%.



Functional Handles



System Mode



**FRONT
INTERFACE**

**Gentle Mode
Sport Mode
Professional Mode**

Before & After



User Experience



Specification

Magnetic Wave :	0-7 Tesla
Total Power :	3000W
RF Power :	120W
Frequency :	F1:1-10Hz, F2:1-100Hz
Pulse Width :	300us
Mode :	Model-I (Sport mode) Model-II (Gentle mode) Model-III (Professional mode)
Screen :	10.4 inch
Treatment Handle :	I-B1, II-B2
Machine Size :	1200mm*420mm*550mm
Package Size :	1230mm*580mm*820mm
Net Weight :	63kg
Gross Weight :	97kg



+91-94460-01470

Business

+91-94460-01472

Support

+91-94460-01470

WhatsApp

info@evolutionmed.com

Email

www.evolutionmed.com

Website