



# HIFEM Chair

Rebuild Strength,  
Reshape Confidence!



30 Min  
Treatment  
Equals  
20,000  
Workout

# HIFEM Chair

## New Technology for Improving Urinary Incontinence

A breakthrough treatment for incontinence and confidence. This unique technology revolutionizes women's intimate health and wellness by providing those suffering from incontinence with a completely non-invasive option.

### How does it work?

Utilize high intensity focused electromagnetic techniques to stimulate the deep muscles of the pelvic floor and restore neuromuscular control. A single focused electromagnetic wave energy brings the greatest contraction on thousands of pelvic floor muscles, which is a very important muscle for incontinent patients. It makes the pelvic floor muscles perform hyper-contraction exercises (20,000 pelvic floor muscle contractions every 30 minutes).



- It use HIEMT (High-Intensity Focused Electromagnetic Wave) technology can target motor nerve cells to stimulate muscle contraction with high intensity, which deeply reshape the internal structure of muscles, strengthen the muscles, increase muscle density and volume.
- Muscle extreme contraction requires a lot of energy supply, then the fat cells next to the muscle will be consumed, so it can help reduce the fat effectively.
- It makes the body lines more prominent, "Mermaid Line" and "The vest Line" clearer.
- It can be said that "you can do exercise while lying down." In 30mins treatment, it equals 20000 workouts.




### Two treatment handles for targeted body parts training



**BELLY, LEGS  
ARMS**



**BELLY, BUTTOCKS  
THIGHS**



## Tired of the Constant Urge to Go?

Pelvic floor muscles are a layer of muscles that support the pelvic organs. Strong pelvic floor muscles provide control over the bladder and bowel. When these muscles are weakened, they do not fully support the internal organs, which may lead to difficulty controlling the release of urine.

## Are you suffering from these afflictions?



Dropped Well-being  
Dissonance in Couple's Life




Leaking urine while laughing, Coughing or Sneezing




Decreased Pelvic Floor Function  
Early Aging  
Menopause



## Traditional Exercise: Kegel Exercises & Pelvic Exercise



1. Preventing urinary incontinence and vaginal prolapse
2. Preventing uterine prolapse
3. Strengthening the pubococcygeal muscles



1. Male Prostate Pain
2. Benign prostatic hyperplasia
3. Prostatitis Magna Erectile Dysfunction

## Traditional Exercise V/S HIFEM Chair Treatment

90% of people cannot exercise deep muscles properly

Unable to know whether the pelvic floor muscles are contracting correctly

The contraction frequency, intensity and amplitude are limited

1

100% contraction of deep pelvic floor muscles.


2

Magnetic energy accurately focuses the pelvic floor muscle penetration depth up to 10 cm.

3

Three-stage treatment mode is effective for training and strengthening.

The strength and frequency of pelvic floor muscle contraction is limited



Ultra-limit contraction of pelvic floor muscles and maintain for a few seconds

- Improve urine leakage and frequent urination.
- One course of treatment = 20,000 Kegel exercises
- Health care of the pelvic cavity
- Sex life is no longer affected
- No need to take off pants, expensive invasive treatment
- Both men and women can use





## Android Smarter & Easier To Operate

- Beauticians Administration
- Patients Administration
- Client Information Storage
- Treatment Date Memory System
- Rental Model



## FEATURES

01 Improve Sexual Sensitivity & Function

02 Prevention and Improvement Stress Urinary Incontinence

03 Postpartum recovery improve vaginal health

04 Strengthen the pubic tail skeletal muscle (PC muscle)

05 Increase blood in the pelvic floor flow & nerve function

06 Increase the vagina after menopause wall thickness & lubricity

95% of patients have reported significant improvement in the quality of life.

### Non-invasive

HIFEM Chair treatment is non-invasive and do not requires and needles at all. HIFEM Energy stimulates the entire pelvic floor area and helps with restoring neuromuscular control.

### Effective

100% contraction of deep pelvic floor muscles, magnetic energy accurately focuses on the pelvic floor muscles with a penetration depth of 10cm.

### Comfortable

HIFEM Chair treatment is comfortable. You will experience a little tingling sensation throughout the 30 minutes at your pelvic area. No numbing cream is required.

## Specification

Magnetic Wave (Energy) : 0-7 Tesla | Voltage : 220V/110V | Output Power : 3200W  
 Frequency : F1:1-10Hz F2:1-100Hz | Pulse Width : 300us | Mode: Sport, Gentle, Professional  
 Screen : 10.4 inch | Machine Size : 830 x 715 x 780 mm



**evolution-med**  
service with a difference



+91-94460-01470  
+91-94460-01472



info@evolutionmed.com  
www.evolutionmed.com