

www.evolutionmed.com

3 IN 1 SYSTEM 8D HIFU

High-Intensity Focused Ultrasound
Skin Lifting System with
Three Treatment Handles.



Treatments

Multi Line: **Body**

Single Line: **Face, Neck, Jaw**

Single Dot: **Around Eyes & Mouth**




8D HIFU

3 IN 1
SYSTEM

High-Intensity Focused Ultrasound Skin Lifting System

8D HIFU is Non Surgical Face Lift Treatment- 8D HIFU (High Intensity Focused Ultrasound) Technology, the aging effects will show on the skin like the wrinkles and the eye-bags. It is the signals to alert you your collagen and elastin are losing from your body. The collagen is the key structural component to support your skin.



8D HIFU is the new technical for face lifting and tightening. 8D HIFU is delivered into the deep skin and cause thermal coagulation to SMAS (Superficial Muscular Aponeurotic System). It can provide the skin tightening and rejuvenation happens gradually with the time. 8D HIFU technology device gives effective, safe and painless skin tightening.

Specifications

Frequency : 4MHz
Power : 350W
Treatment Depth : DS1.5, DS3.0, DS4.5mm (Standard)
DS6.0, DS8.0, DS10.0, DS13.0mm,
DS16.0mm (Optional)
Treatment Length : 8-25mm
Width : 0--10mm
Energy : 0.1-- 2.0J
Pitch : 1.0mm-2.5mm
Clearance : 1.0mm-2.0mm
Power Source : AC 110V~220V--50Hz/60Hz



+91-94460-01470
+91-94460-01472



info@evolutionmed.com
www.evolutionmed.com

Kerala | Tamil Nadu | Maharashtra | Dubai

Advantages

- Completely non-invasive and non- ablative to the epidermis layer.
- No down time, and make- up can be wearied only in 15 min after the treatment.
- Fast operation time , full face and neck only takes 40-60 minutes.
- Long lasting result, 2-3 year can last only with one treatment.
- Flexible treatment, one single treatment achieves instant result. High quality, low consumables.

Treatment Effect

- Strengthen and tighten relaxed skin with remarkable tightening effect.
- Only one treatment is needed to smooth wrinkles and improve them comprehensively.
- Make skin texture more delicate and glossy.
- Filling in shrinkage and loss of collagen vacancies to enhance the relaxation profile.
- Remove the traces of age such as roughness, and Skin rejuvenation.
- Make skin tender, white and youthful.

Before & After

